

A REAL STEP FORWARD IN BIOSECURITY

HYSOLV
FOOTCHECK
THE PROFESSIONALS CHOICE



- **Simple, Economical and Easy-To-Use**
- **Accurate and correct concentration every time**
- **Audit tracking for recording and responsible use**



A REAL STEP FORWARD IN BIOSECURITY

HYSOLV FOOTCHECK THE PROFESSIONALS CHOICE



Measuring chambers

Disinfectant poured direct from container, ensures correct dilution, minimises waste

0.5%
chamber

1%
chamber

2%
chamber

Possible concentrations

0.5%, 1.0%, 2.0%, 2.5%,
3.0% & 4.0%*


* Fill level B required

Durable HDPE

Manufactured in durable HDPE, with UV protective additive, designed to be robust and easily cleaned

Foot sole rub bars

Foot safety ribbing.
Aid to help avoid slipping




Audit stickers

Panel for peelable audit stickers for simple recording of contents and frequency of changes

Calibrated fill level (15 litre)

Ensures precise dilution avoiding excessive use of expensive concentrate



Ready-to-use record sheet -
Download from the website

Removable hinged lid

To prevent contamination, evaporation or dilution by rain